



# PC-1000/Laser 1000 Panoramic/Cephalometric X-ray Machine



Performance, versatility and value make the PC-1000/Laser 1000 #1 in the nation

Choosing the panoramic/cephalometric x-ray machine that's right for your practice is an important decision. Key considerations include enhancing patient care, increasing profitability and keeping pace with technology.



The Panoramic PC-1000/Laser 1000 is the choice of more dental practitioners than any other panoramic/cephalometric x-ray machine in the nation. It has also received the prestigious *DentalTown Magazine* Townie Choice Award for best panoramic/cephalometric x-ray multiple years in a row.

Recognizing that each dental practice is different, and that circumstances can change as your practice grows, Panoramic designed the PC-1000/Laser 1000 to provide exceptional quality and value today, as well as far into the future.

## PC-1000/Laser 1000 Highlights:

- Consistent, high-quality panoramic/cephalometric x-ray images, exceeding those of machines costing far more.
- The market's simplest control panel and overall ease-of-use.
- Legendary durability backed by the industry's longest and most inclusive warranty.
- Remarkable value, which includes shipping, installation and staff training.
- Responsive customer service and expert technical support.

## Learn more about the PC-1000/Laser 1000 today.

See reverse side for product details and call us at 1-800-654-2027 or visit our web site at [www.pancorp.com](http://www.pancorp.com). Then make the decision that's right for your practice.



# PC-1000/Laser 1000 Panoramic/Cephalometric X-ray Machine

*Performance, versatility and value set the PC-1000/Laser 1000 apart from all other panoramic/cephalometric x-ray machines, making it the best-selling model on the market.*

## Performance

- Consistent, high-quality images.
  - Our 0.5mm x 0.5mm **focal spot** is among the smallest in the industry, producing a clearer, sharper image.
  - Our 12mm anterior and 41mm posterior **focal trough** is among the largest in the industry, making patient positioning easier and more repeatable to ensure accuracy.
- Simple operation.
  - The PC-1000/Laser 1000 features the **simplest control panel and overall ease-of-use** on the market, reducing the risk of error and minimizing learning curves.
  - The **fixed bite guide** eliminates guesswork and the easy-to-use adjustment guides help ensure proper patient positioning.
  - Integrated **soft tissue filter**.
- Quality workmanship.
  - **Individually manufactured** at our Midwest headquarters for lasting reliability.
- Exceptional service and technical support.
  - **A nationwide network of trained field service technicians** to ensure your needs are promptly met.
  - **Free telephone and Internet support** from knowledgeable factory service technicians.
- Legendary durability backed by **the industry's longest and most inclusive warranty**.

## Versatility

- **Compatible with virtually any phosphor plate** package on the market.
- **Freestanding** with no wall or floor mounting required, allowing ease of installation, maintenance and transportability.
- **Control panel** can be mounted on either right or left side of unit, adapting easily to room layout.
- PC-1000/Laser 1000 **eliminates the need for two separate machines**.
- The integrated Cephalostat is **adjustable in 45° increments**, which facilitates seven common position views:
  - Lateral Skull
  - Lateral – Oblique
  - Posterior – Anterior
  - Anterior – Posterior
  - Occipitontal (Water's View)
  - Reverse – Towne's
  - Submentovertex
- The PC-1000/Laser 1000 provides the **variety of angles, views and measurements** you want.

## Value

- **Overall affordability** with performance topping far more costly units.
- **Flexible financing, purchase and leasing options**, including our innovative "Income-Per-Image" program.



**"One Plan. One Pan." – Flexibility and Assurance.** Anyone who has invested in a new technology only to learn all too soon that a newer product does it faster, better and cheaper, will appreciate Panoramic's "One Plan. One Pan." commitment. Recognizing that each dental practice is different and evolves over time, **"One Plan. One Pan." guides our products and services to anticipate and adapt to your clinical, technological and financial needs.** It empowers you to adopt new technologies at your own pace – quickly, affordably and confidently. "One Plan. One Pan." embodies the founding principles that have made Panoramic the market leader: listen to the dentist and deliver the simplest solutions at the best value.



Panoramic Corporation

Contact one of our Product Specialists to help you evaluate your needs and customize a solution.

phone 800.654.2027 www.pancorp.com sales@pancorp.com