



PC-1000 Panoramic X-ray Machine

Performance, versatility and value make the PC-1000 #1 in the nation

Choosing the panoramic x-ray machine that's right for your practice is an important decision. Key considerations include enhancing patient care, increasing profitability and keeping pace with technology.



The Panoramic PC-1000 is the choice of more dental practitioners than any other panoramic x-ray machine in the nation. It has also received the prestigious *DentalTown Magazine Townie Choice Award* for best panoramic x-ray multiple years in a row.

Recognizing that each dental practice is different, and that circumstances can change as your practice grows, Panoramic designed the PC-1000 to provide exceptional quality and value today, as well as far into the future.

PC-1000 Highlights:

- Consistent, high-quality panoramic x-ray images, exceeding those of machines costing far more.
- Easy, affordable conversion to state-of-the-art digital technology.
- The market's simplest control panel and overall ease-of-use.
- Legendary durability backed by the industry's longest and most inclusive warranty.
- Remarkable value, which includes shipping, installation and staff training.
- Responsive customer service and expert technical support.

Learn more about the PC-1000 today.

See reverse side for product details and call us at 1-800-654-2027 or visit our web site at www.pancorp.com or www.pancorpdigital.com. Then make the decision that's right for your practice.



PC-1000 Panoramic X-ray Machine

Performance, versatility and value set the PC-1000 apart from all other panoramic x-ray machines, making it the best-selling model on the market.

Performance

- Consistent, high-quality images.
 - Our 0.5mm x 0.5mm **focal spot** is among the smallest in the industry, producing a clearer, sharper image.
 - Our 12mm anterior and 41mm posterior **focal trough** is among the largest in the industry, making patient positioning easier and more repeatable to ensure accuracy.
- Simple operation.
 - The PC-1000 features the **simplest control panel and overall ease-of-use** on the market, reducing the risk of error and minimizing learning curves.
 - The **fixed bite guide** eliminates guesswork and the easy-to-use adjustment guides help ensure proper patient positioning.
- Quality workmanship.
 - **Individually manufactured** at our Midwest headquarters for lasting reliability.
- Exceptional service and technical support.
 - **A nationwide network of trained field service technicians** to ensure your needs are promptly met.
 - **Free telephone and Internet support** from knowledgeable factory service technicians.
- Legendary durability backed by **the industry's longest and most inclusive warranty.**

Versatility

- Digital-ready, enabling you to adopt new technologies at your own pace.
 - **Becomes a state-of-the-art direct-digital system** with the 1000-DR Digital Conversion Package.
 - **Compatible with virtually any phosphor plate** package on the market.
- **Freestanding** with no wall or floor mounting required, allowing ease of installation, maintenance and transportability.
- **Control panel** can be mounted on either right or left side of unit, adapting easily to room layout.
- Available with a **cephalometric unit** (also compatible with phosphor plate packages), eliminating need for two separate machines.

Value


- **Overall affordability** with performance topping far more costly units.
- **Flexible financing, purchase and leasing options**, including our innovative "Income-Per-Image" program.

Cephalometric Capabilities Available

The PC-1000/Laser 1000 Cephalometric attachment adds an important capability to your practice.



"One Plan. One Pan." – Flexibility and Assurance. Anyone who has invested in a new technology only to learn all too soon that a newer product does it faster, better and cheaper, will appreciate Panoramic's "One Plan. One Pan." commitment. Recognizing that each dental practice is different and evolves over time, **"One Plan. One Pan." guides our products and services to anticipate and adapt to your clinical, technological and financial needs.** It empowers you to adopt new technologies at your own pace – quickly, affordably and confidently. "One Plan. One Pan." embodies the founding principles that have made Panoramic the market leader: listen to the dentist and deliver the simplest solutions at the best value.



Panoramic Corporation

Contact one of our Product Specialists to help you evaluate your needs and customize a solution.
phone 800.654.2027 www.pancorp.com www.pancorpdigital.com sales@pancorp.com