


## PC-4000



### ■ Specifications

Completely Free Standing

Physical Dimensions: 42" x 35"

Recommended Minimum  
Working Space: 48" x 48"

Focal Spot: 0.5mm x 0.5mm

Focal Trough:  
Anterior Region - 12mm  
Posterior Region - 41mm

Rated kVp: 70-90

Tube Current: 6 mA

Power Requirements:  
115 V, 20 Amp Dedicated Service

Exposure Time: 12 Seconds

Image Size Captured:  
6" x 12" (PAN and TMJ)

Weight: 415 lbs.

Maximum Height: 91"

### ■ Space Requirements

Physical Dimensions: 42" x 35"

Recommended Working Space: 48" x 48"

### ■ Electrical Requirements

115 Volt Grounded Outlet

20 Amp Service